MOBILITY UK
Finding Solutions

LifeGlider®

LifeGlider helps you walk safely.

If you or a loved one suffers with mobility issues, then you already know the difficulties faced and the upset for everyone if there is a fall.



Mobility.org.uk



Other walking aids like frames, sticks, and crutches can increase the risk of falling, leading to injuries, harm and distress.

LifeGlider is Different.

LifeGlider helps people who can walk a little but struggle with balance, strength, stability or orientation.

It supports many conditions including Stroke, MS, Parkinsons, FND, CES, Cerebral Palsy, Ataxia, Neuropathy and age related mobility loss.

It keeps you upright, secure and able to walk with confidence.

- Walk Safely
- Natural Standing Position
- Better Posture
- Great For Rehabilitation
- Perching Seat to Rest



LifeGlider®

MOBILITY UK
Finding Solutions

LifeGlider®

LifeGlider is the first fall safe hands-free upright mobility device.

*** On flat safe surface, you can also use the parallel hand grips for support.**

1. Secures centre of gravity so users won't fall when they lose their balance or lose their grip.
2. Hands optional use frees hands to do other activiters, cook, hold devices, dance.
3. Allows eye level social contact, look ahead vs looking at the ground.
4. Enables standing and walking for longer periods to strengthen aiding rehabilitation.
5. Slows or prevents the need for a wheelchair and other mobility devices.
6. Gives a great mental health boost and dramatically improves psychological well being.



For more information or to book a video demo.



National Helpline
0300 555 3555

WhatsApp Chat
07777 300 769

Or send us an email
lg@mobility.org.uk



Order Yours Today
Mobility.org.uk

LifeGlider®

MOBILITY UK
Finding Solutions

Is the LifeGlider right for you or your patient?



RECOMMENDED

Can walk a little without or with assistance.
Balance issues due to age, injury, or neurological conditions.
e.g. Cerebral Palsy, Multiple Scleroris Parkinsons, FND, CES, Stroke

NOT RECOMMENDED

- Advance Dementia.
- Unable to Walk or Stand.
- Visually Impaired
- Has Hemiplegia
- Extreme Spasms
-


IMPORTANT CRITERIA FOR USE

User must be able to stand for one minute, long enough to secure themselves into the LifeGlider

Must have sufficient grip strength in at least one hand to apply the brake.

The secret to fall prevention is to secure the center of gravity

Traditional walkers shift the center of gravity, increasing the risk of falling.



1. Located at the small of the back.
2. Ensures center of gravity remains at the core.
3. Improves posture, gait and balance.

LifeGlider®

MOBILITY UK
Finding Solutions

Recommended by Experts



I regret not finding the LifeGlider sooner as it would have been absolutely amazing for my mum following her stroke. It has been great for my daughter who has CES and FND. I think it is revolutionary.

Martyn, founder of Mobility UK

mobility.org.uk | 0300 555 3555

MOBILITY UK
Finding Solutions



I help people to train in a range of equipment and the LifeGlider is an intrinsic part of our routines.

Matt has been offering and campaigning for inclusive exercise for all and recently opened the Angel Centre in Colchester, a dedicated inclusive gym.

Matt Brinkley | therapygyms.co.uk

Therapy Gyms
Making fitness accessible



So happy to have found this and provided it to our patient who has gained new freedom. We offer physio in the LifeGlider and love the results.

Bespoke Rehabilitation Physiotherapy.
Maximising outcomes: needs & quality driven.

Lisa - reachrehab.co.uk

Reach®
Rehab Physiotherapy



When I discovered the LifeGlider, I knew instantly that it would help my Mum and many others.


Ewan is a leading advocate for regaining freedom and independence. He offers expert advice on the LifeGlider, group demonstrations and is well known across the industry.

Ewan Phillips - Regain Mobility

Regain
MOBILITY

Customer Stories

Walking handsfree, safely, changes lives.




McKenna Baker

Cerebral Palsy

I have had CP since I was born.

With LifeGlider, walking is easier and I am even dancing.

I feel safe and can even walk my daughter to school and dance with my Mum.




Kelly Gregory

Friedreich Ataxia

I've been walking with Ozzie (yes I name my LG) and getting stronger every day.

Watching back the videos we make shows me my progress and I am walking more upright and gaining better balance.

Insta: @conqueringfa




Lana Kurasidze

Stage 4 Brain Tumor

At 21, 15 years ago, doctors told me I only had three months to live. I lost the ability to walk and have been fighting back ever since.

I went from Wheelchair to Walking. I even went shopping for the first time in 11 years.

It's giving me my life back.









A Neurologist's View

Dr. Adornato, a neurologist with over 45 years of experience, praises the LifeGlider for its exceptional safety and stability for patients with impaired mobility. Unlike conventional walkers, the LifeGlider prevents falls, allows hands-free use, and supports various activities.

It not only enhances mobility but also aids in rehabilitation, improving users' strength, confidence, and quality of life.

Bruce Adornato MD
Adjunct Clinical Professor of Neurology and Neurological Sciences, Stanford School of Medicine
https://med.stanford.edu/profiles/bruce-adornato

Stanford MEDICINE



Width
21 inches
53 cm

1. Adjustable belt to secure your center of gravity.
2. Brake designed to give you greater control can be installed on either side.
3. Folds for easy transport or storage.
4. Durable high-strength steel construction.
5. Seat is angled precisely to shift weight off your feet when needed.
6. Padded hand grip; use lightly when walking, if you like, to feel extra secure.
7. Adjustable height for maximum comfort.
8. Specially-designed wheels support weight rating.

Available in Three Sizes

Sizes	SMALL	STANDARD	TALL
Height	1.37-1.63m	1.5-1.73m	1.73-1.93m
Inseam	61-68 cm	71-76 cm	76-96 cm
Max Weight	90.7kg	124.7kg	124.7kg

* Maximum weight support includes the weight of objects you are carrying.

